

bepalingen	matrix	scope	methode	externe kwal. controle
aceton	serum	KF.TOX.01	GC-FID	SKML + LGC
amitriptyline	serum	KF.TDM.01	LC-MS/MS	LGC
carbamazepine	serum	KF.TDM.02	LC-MS/MS	LGC
ciclosporine	EDTA volbloed	KF.TDM.02	LC-MS/MS	LGC
clomipramine	serum	KF.TDM.01	LC-MS/MS	LGC
clozapine	serum	KF.TDM.01	LC-MS/MS	LGC
coffeine	serum	KF.TDM.01	HPLC-DAD	SKML
desipramine	serum	KF.TDM.01	LC-MS/MS	LGC
desmethylclomipramine	serum	KF.TDM.01	LC-MS/MS	LGC
digoxine	serum	KF.TDM.03 (advies)		LGC
DOA	urine	KF.TOX.03	EIA	SKML
ethanol	serum	KF.TOX.01	GC-FID	SKML + LGC
ethanol	serum	KF.TDM.03 (advies)		SKML + LGC
ethyleenglycol	serum	KF.TOX.01	GC-FID	SKML + LGC
everolimus	EDTA volbloed	KF.TDM.02	LC-MS/MS	LGC
fenobarbital	serum	KF.TDM.03 (advies)		LGC
fenytoine	serum	KF.TDM.03 (advies)		LGC
gentamicine	serum	KF.TDM.03 (advies)		SKML + LGC
GHB	serum	KF.TOX.01	GC-FID	SKML + LGC
imipramine	serum	KF.TDM.01	LC-MS/MS	LGC
lamotrigine	serum	KF.TDM.02	LC-MS/MS	LGC
lithium	serum	KF.TDM.03 (advies)		SKML
levetiracetam	serum	KF.TDM.02	LC-MS/MS	LGC
methanol	serum	KF.TOX.01	GC-FID	SKML + LGC
norclozapine	serum	KF.TDM.01	LC-MS/MS	LGC
nortriptyline	serum	KF.TDM.01	LC-MS/MS	LGC
paracetamol	serum	KF.TDM.03 (advies)		SKML + LGC
sirolimus	EDTA volbloed	KF.TDM.02	LC-MS/MS	LGC
STIP	serum	KF.TOX.01	HPLC-DAD	SKML
tacrolimus	EDTA volbloed	KF.TDM.02	LC-MS/MS	LGC
theofylline	serum	KF.TDM.03 (advies)		LGC
thiocyanaat	serum	KF.TOX.01	UV-VIS	spike
topiramaat	serum	KF.TDM.02	LC-MS/MS	LGC
valproïnezuur	serum	KF.TDM.03 (advies)		LGC
vancomycine	serum	KF.TDM.03 (advies)		SKML + LGC

cortisol	urine/speeksel	CH.KCH.07	LC-MS/MS	SKML
cortison	urine/speeksel	CH.KCH.07	LC-MS/MS	SKML
HbA1c + Hb variant	EDTA volbloed	CH.KCH.06	HPLC	SKML
metanefrines	EDTA plasma	CH.KCH.07	LC-MS/MS	SKML
methylmalonzuur	serum	CH.KCH.07	LC-MS/MS	SKML
DHEAS (steroiden)	serum	CH.KCH.07	LC-MS/MS	SKML + Instand
androsteneendion (steroiden)	serum	CH.KCH.07	LC-MS/MS	SKML + Instand
11-deoxycortizol (steroiden)	serum	CH.KCH.07	LC-MS/MS	SKML + Instand
17-hydroxy progesteron (steroiden)	serum	CH.KCH.07	LC-MS/MS	SKML + Instand
vitamine B1	EDTA volbloed	CH.KCH.07	LC-MS/MS	SKML
vitamine B6	EDTA volbloed	CH.KCH.07	LC-MS/MS	SKML
Vitamine A	serum	CH.KCH.07	LC-MS/MS	SKML
Vitamine E	serum	CH.KCH.07	LC-MS/MS	SKML